Aesculap Spine

- (GB)** Instructions for use/Technical description
activO frame system
- (USA)** **(D)** Gebrauchsanweisung/Technische Beschreibung
activO-Rahmensystem
- (F)** Mode d'emploi/Description technique
Système de cadre activO
- (E)** Instrucciones de manejo/Descripción técnica
activO sistema separador con marco
- (I)** Istruzioni per l'uso/Descrizione tecnica
Sistema telaio activO
- (P)** Instruções de utilização/Descrição técnica
Sistema de armação activO
- (NL)** Gebruiksaanwijzing/Technische beschrijving
activO-Framesysteem

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
activO frame system

Legend

- A Connecting piece with quick-release coupling
- B Coupling bar
- C Coupling piece short
- D Radiolucent valve
- E Fixation screws
- F Screwdriver for fixation screws
- G Fastening clamp
- H T-wrench
- I OP field illumination bar
- J UNITRAC universal holder
- K Forceps for changing the blades

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bb Braun.com

Intended use

The Aesculap activO frame system is used to keep free a ventral or ventrolateral approach to the spine. With this system, valves can be placed at any position on the ring and fastened on the spine, using fixation screws.

Illumination and optics elements can be adapted, too. The frame can be fastened to the operating table (see instructions for use TA011611 or TA009150 respectively). In the latter configuration, fixation with fixation screws on the spine is not required.

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Available sizes

Art. no.	Designation
BV399R	Forceps for changing the blades
BV858R	Connecting piece with quick-release coupling
BV859R	Coupling piece short
BV860R	Coupling bar 150 mm
BV861R	Coupling bar 200 mm
BV862P	Radiolucent valve: Size S-green
BV863P	Radiolucent valve: Size S-yellow
BV864P	Radiolucent valve: Size L-blue
BV865SU	Fixation screw: Size S-green
BV866SU	Fixation screw: Size M-yellow
BV867SU	Fixation screw: Size L-blue
BV868R	Screwdriver for fixation screws
BV869R	Fastening clamp
BV899R	T-wrench
OP803	OP field illumination bar 100 mm

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

Note

Use of the activO frame system requires precise knowledge about the biomechanical situation at the spine.

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ Do not use products that are in open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.

The products BV865SU, BV866SU and BV867SU are sterilized by radiation and are supplied in sterile packaging.

The product must not be reused.

- ▶ Do not re-sterilize the product.

Aesculap®

activO frame system

Safe operation

Mounting the radiolucent valves

Injuries to nerve roots and spinal cord due to mismatch between fixation screw and valve sizes!

- Use the color-coded fixation screws matching the color of the valve.
- Screw in the fixation screw under X-ray monitoring.

Injury to soft tissue parts or vessels, caused by fixation screws!

- Position the valve without a fixation screw. Only when the valve is correctly positioned, insert and tighten the fixation screw.

Breakage, deformation or loosening of the fixation screw due to manipulating the valve after insertion of the fixation screw, resulting in complication or prevention of fixation screw removal from the valve!

- Place the valve in the correct position before inserting the fixation screw.
- When using a valve with a fixation screw: Affix the valve to the associated frame using the coupling bar.

Valve breakage!

- Avoid knocks on the valve or on components connected to the valve.
- Do not impel the fixation screws with a mallet under any circumstances.

- ▶ Select a valve **D** of appropriate length (S, M or L) for the patients.
- ▶ Valve **D** can be fixed to the vertebral body with a fitting fixation screw **E**. Be careful to use a valve and screw of the same size (see color code).
- ▶ Apply screwdriver **F** on the hexagon head of fixation screw **E**.
- ▶ Insert and tighten fixation screw **E** in valve **D**, see Fig. 1.


Fig. 1

Loosening the fixation screw

- ▶ Apply screwdriver **F** on fixation screw **E** and unscrew counterclockwise.
- ▶ Remove screwdriver **F** with fixation screw **E** from the operating field.

Note

The fixation screw **E** is fixed on the screwdriver **F**.

Adapting the coupling bars

Damage to soft tissue parts!
Loss of approach due to incorrect tightening or mounting of the coupling bars!

- ▶ Tighten the nut of the coupling bar, then verify that the angular position cannot be moved anymore.
- ▶ Always engage the valve coupling until it clicks.

- ▶ Depending on application, mount coupling piece **C** or coupling bar **B** on radiolucent valve **D**.

Note

Optionally, coupling bars **B** can be mounted at an offset angle of 90°.

Aesculap®

activO frame system

- If necessary mount universal holder **J** on the respective coupling.
- To extent the construct, install connecting piece with quick-release coupling **A** on the respective coupling, see Fig. 2.


Fig. 2

Adjusting the inclination angle of the coupling bars

Note

*Coupling bars **B** can be adjusted to an inclination of between 15° upwards and 90° downwards.*

- To loosen the coupling lock, apply T-wrench **H** on hexagon **1** and untighten counterclockwise.

- Move T-wrench **H** to the required angle and tighten hexagon **1** turning the tool clockwise, see Fig. 3.
Coupling bar **B** is locked again.


Fig. 3

Releasing the quick-release coupling

- Grip the quick-release coupling **1** of the connecting piece **A** or the UNITRAC universal holder **J** on the handle part and pull backwards, see Fig. 4.


Fig. 4

Note

*Optionally, optics holder **RT081R** can be fixed to UNITRAC universal holder **J** and connecting piece with quick-release coupling **A**.*

- Fasten connecting piece with quick-release coupling **A** on the frame, using fastening clamp **G**, and install optics holder RT081R (see instructions for use TA009907).

Uncoupling the radiolucent valve

Note

Radiolucent valve **D** can be uncoupled in two different ways.

- Grasp the coupling and press with the thumb on pin **1** of valve **D**.
– or –
- Push forceps **K** with the jaw side **2** under the coupling and use the forceps bolt **3** to push them onto the pins **1** of the valve **D**, see Fig. 5.


Fig. 5

Mounting the fastening clamp on the frame

WARNING

Loss of approach due to incorrect mounting and closing of the fastening clamp!

- Firmly tighten the fastening nut.
- Always engage the coupling bars and the frame until they click.

CAUTION

Loss of the fastening nut and the upper clamping element caused by complete removal of the fastening nut!

- Only loosen the fastening nut by 2–3 turns.

- Loosen fastening nut **1** by 2–3 counterclockwise turns.
- Click the fastening clip **G** into the desired position onto the frame **4** (BV814R/BV815R), see Fig. 6. Fastening clamp **G** is fixed on frame **4**.


Fig. 6

Aesculap®

activO frame system

- Click the coupling bar **B** or connecting piece with the quick-release coupling **A** into the jaw piece of the fastening clamp **G** and place it in the desired position, see Fig. 7.


Fig. 7

- Tighten fastening nut **1** until the clamp is fully closed. If necessary tighten further with T-wrench **H**, see Fig. 8.


Fig. 8

Dismounting the fastening clamp from the frame

- Loosen fastening nut **1** by 2–3 counterclockwise turns, using T-wrench **H** if necessary.
- Unclick coupling bar **B** or connecting piece with quick-release coupling **A** from the jaw part of fastening clamp **G**.
- Slide off fastening clamp **G** from the frame by applying lateral pressure on the clamp, see Fig. 9.


Fig. 9

Mounting the OP field illumination bar

Burns from the hot end of optical fiber!

- ▶ When using the product with a light source, make certain that the end of the optical cable does not touch human tissue or highly inflammable materials. Keep a distance of at least 10 mm.
- ▶ Only use light sources of a power rating of up to 180 W.

- ▶ Affix the surgical lighting rod **I** directly with the fastening clip **G** or, using the bracket BV020R and the bracket for OP surgical lighting FG410R affix onto the coupling bars **A**, **B** or **C**.

The OP field illumination bar has three different light cable connectors, see Fig. 10.

- ▶ Unscrew connectors that are not required.


Fig. 10

Disassembling/Assembling the fastening clamp

Disassembling

- ▶ Remove fastening nut **1** from threaded pin **3** by turning it counterclockwise.
- ▶ Remove jaw part **2**, see Fig. 11.

Fig. 11

Assembling

- ▶ Put jaw part **2** with its hole over threaded pin **3**.
- ▶ Screw on fastening nut **1** clockwise, with the hexagon facing upwards.

Aesculap®

activO frame system

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Single-use products

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► Do not reprocess the product!

Note

Products for single use are BV865SU, BV866SU and BV867SU.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the instructions for use (BV869R).

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Risk to patients!

- BV858R, BV860R, BV861R, RT073R, RT081R, FG410R and BV020R should only be cleaned mechanically!

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable temperature of 55 °C.

Aesculap®

activO frame system

- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ For BV399R, BV842R, BV862P, BV863P, BV864P, BV868R, BV859R, BV869R, BV899R, FK166R, RT092R, RT069R, and OP803	■ Suitable cleaning brush:, e.g. TA011944 ■ 20 ml disposable syringe ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ For BV814R/BV815R and BV843R	■ Suitable cleaning brush:, e.g. TA011944 ■ 20 ml disposable syringe ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ For BV399R, BV859R and FK166R	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Mechanical neutral or mild alkaline cleaning and thermal disinfecting ■ For OP803	■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Aesculap®

activO frame system

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ For BV862P, BV863P, BV864P, BV868R, BV869R, BV842R, RT069R and BV899RP	■ Suitable cleaning brush; e.g. TA011944 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical neutral or mild alkaline cleaning and thermal disinfection ■ For RT092R	■ Suitable cleaning brush; e.g. TA011944 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical neutral or mild alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ For BV814R/BV815R, BV843R, BV858R, BV860R, BV861R, RT073R, RT081R, FT410R and BV020R	■ Suitable cleaning brush; e.g. TA011944 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.

► Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

Aesculap®

activO frame system

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

Aesculap®

activO frame system

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

activO frame system

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: ■ Concentrate: - pH neutral - <5 % anionic surfactant ■ 0.5 %* working solution Mildly alkaline: ■ Concentrate: - pH = 9.5 - <5 % anionic surfactant ■ 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Aesculap®

activO frame system

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

activO frame system

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: ■ Concentrate: - pH neutral - <5 % anionic surfactant ■ 0.5 %* working solution Mildly alkaline: ■ Concentrate: - pH = 9.5 - <5 % anionic surfactant ■ 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

The product can be sterilized either in disassembled or in assembled condition.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min or as per locally approved policies and procedures
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Aesculap®

activO frame system

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/ 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042

Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany